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Index

Indexation to the GvOC, why ?	p.2
The GPG Clamp	3
Method of use of the GPG	3
Placement of the retractors	4
Guide wire: important notice	5
Introducing the guide wire: caution	5

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
Indexation to the GvOC, why ?

Incorrect placement of a glenoidal implant directly influences the clinical and radiographic results of total shoulder arthroplasty (1).

The placement of a glenoidal implant for total shoulder prosthesis is made difficult (a) by a superficial aspect of the glenoid bone which is not always reliable to determine the deep geometry of the scapular neck, and (b) by the erosion of the arthritic glenoid which alters the interpretation of its native anatomical orientation.

Several anatomical studies have shown that the plane of the scapular blade is not a reliable landmark for determining the retroversion, tilt, rotation and offset planes of the pre-pathological native glenoid. Especially since this reference plane evolves with the age of the patient (2). The authors show that the plane of the Glenoid vault Outer Cortex (Gvoc) is comparatively more reliable, and especially does not seem to evolve with the aging of the patient.

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The GPG clamp *La pince GPG*

Placing the GPG (Glenoid Positioning Guide) clamp on the shoulder blade neck is simple and reproducible.

The GPG clamp allows the installation of a guide wire simultaneously at the centre of the surface of the glenoid bone and in the orientation of the shoulder blade neck.

The jaws of the clamp are positioned on the anterior and posterior sides of the scapular neck at the GvOC level – approximately 15 to 20mm medial to the glenoid rim.
2 guiding tunnels allow to position the pin guide either centered in the middle of the clamp, or shifted vertically by 7mm.

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
During the total shoulder prosthesis surgery, of anatomical or reverse type, the glenoid preparation time will come once the humeral preparation time has been done. This order is necessary in order to allow comfortable access to the glenoid.

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1- The first step on the glenoid will be to determine the center of the glenoid and its orientation: after anterior and posterior capsulotomy, and labrum resection, trace the vertical and the transversal axis of the glenoid surface with a diathermy in order to visualize its center point (fig.1).

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The positioning of the retractors is of paramount importance to allow the correct placement of the GPG on the neck of the scapula.

- (1) Place a hohmann retractor at the postero-superior part of the glenoid.
- (2) Place another long S retractor on the anterior part of the glenoid below the sub-scapularis.
- (3) Place the Trouilloud retractor at the postero-inferior part of the glenoid, allowing the humerus to be pushed back.

Depending on whether the operator is preparing an anatomical or reverse TSA:

- Anatomic: locate the equator of the glenoid and mark a position line with a diathermy.
- Reverse: place a mark with a diathermy 15mm from the lower margin of the glenoid.


(4) Leaving the retractors in place, open the posterior capsule -about 2cm long- and introduce a periosteal to the posterior part of the glenoid to completely remove the capsular attachments.

(5) Remove the anterior retractor, and open the anterior capsule flush with the glenoid -about 2cm long- and pass a periosteal to the anterior border of the glenoid to completely remove the capsular attachments.

(6) Introduce the branches of the GPG on both sides of the glenoid into the capsular openings.

Adjust the position of the GPG in the vertical plane so that the guide wire points to the center of the mark made with the diathermy (height adjustment, not angulation).

In the case of a Reverse TSA, consider keeping a downward tilt when securing the GPG on the scapular neck.

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Guide wire: important notice

Note importante : broche de guidage

Guide wires (Ref 33-0225-0310), even genuine, must be controlled when preparing the instrumentation set before sterilization, and if in doubt, before being used during surgery in the centering hole of the GPG clamp:

- place each wire on a perfectly flat surface and roll it with a finger.
- the wire must roll without resistance and without coupling.
- if in doubt, remove the wire from the instrumentation set and replace it with a new wire that has been checked using the same method.

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fig.2

2- Introduce the GPG clamp on either side (anterior and posterior) of the glenoid neck (fig.2), 10 to 15mm internal to the glenoid surface, then clamp the GPG on the Glenoid Vault outer Cortex (GvOC)(fig.3).

fig.3

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3- Introduce the guide wire (ref. 33-0225-0310) into the center hole of the clamp. Most of the times the spontaneous positioning of the GPG clamp will lead the guide wire straight towards the center point of the glenoid surface (fig.4).


fig.4

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Due to a very tight fitting of the guide wire in the holes of the GPG allowing high positioning accuracy, the wire must be perfectly aligned with the axis of the hole and no resistance must be felt by the operator when the wire is first introduced manually.

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4- GPG in place, the version is strictly defined by the positioning of the GPG and the surgeon will have to adapt vertically the entry point and its vertical angulation to the requirements of the on-going surgical procedure:

4a: for the correct placement of a UNIC Reverse glenoid baseplate, the entry point should be 14-15mm above the distal edge of the glenoid surface, AND tilted down by 5 to 10°.

4b: for the placement of a UNIC anatomic glenoid component, the entry point is the center drawn at the cross-road of the vertical and transversal axis (point 1 here-above) in neutral orientation.

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5- Drill the guide wire into the glenoid.
The guide wire is graduated by steps of 5mm from 15 to 30mm of introduction.
Introduce the guide wire through the first cortex and until contact with the second cortex and without drilling it (fig.5).

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
6- Visually measure the introduction of the guide wire by counting the remaining markers out of the glenoid:


- No visible mark = introduction of 30mm minimum,
- 1 visible mark = introduction between 25 and 30mm,
- 2 visible marks = introduction between 20 and 25mm,
- 3 visible marks = introduction between 15 and 20mm,
- 4 visible marks = introduction less than 15mm.

In order to secure the placement of the Glenoid implant, the measured value must be higher than that of the implant keel (22mm for a UNIC standard glenoid base, 17.5mm for a UNIC XS glenoid base).

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7- In cases where the guide wire is not positioned satisfactorily, it is possible to use the second aiming hole of the GPG which is offset by 7mm (fig.7). In case the offset hole is used, the GPG clamp can be used in its "conventional" position for an 7mm proximal offset, or can be "reversed" for an 7mm distal offset (fig.8).

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fig.8

8- Continue introduction of the guide wire with the powertool until fixation into the 2nd cortex.

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

fig.9

Once the guide wire is in place, remove the GPG clamp leaving the guide wire in position (fig.9). The remaining steps of the implant's operating technique will be performed in accordance with the implant's surgical technique manual.

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GPG® Glenoid Positioning Guide

(1) Accuracy of Glenoid Component Placement in Total Shoulder Arthroplasty and its Effect on Clinical and Radiological Outcome in a Retrospective, Longitudinal, Monocentric Open Study. T.Gregory, A.Sankey, B.Augereau, E.Vandenbussche, A.Amis, R.Emery, U.Hanssen. Plos one-Oct 2013 - Volume 8 - Issue 10 - e75791.

(2) Does the Orientation of the Glenoid Rim respect the Scapula Blade evolving with Age compared to the Glenoid Vault Outer Cortex Plane? T.Gregory, M.Artuso. SECEC 2019 Abstract Proceedings - JSES Open Access 3 (2019) 233-260.

(3) Is the Glenoid Vault Outer Cortex a more Accurate Reference Plane than Conventional Methods in Shoulder Arthroplasty? T.Gregory, S.Hurst, L.Merlini, U.Hansen, J.Gregory, R.Emery. DOI: <https://doi.org/10.21203/rs.3.rs-93174/v1>