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GB - EVOLUTIS INSTRUMENTS
INSTRUCTIONS FOR USE
VERSION 2025-03

For the attention of the surgeon and operation theatre staff

1. WARNINGS

EVOLUTIS instruments are to be exclusively used for implantation or explantation of EVOLUTIS prosthesis. They are composed of: cannulated (or not) drills / reamers, tray, scalers, screwdrivers, rasps, burs, trial prosthesis, cutting guide, handles, extractors, impactors, grippers, guide pins for knee, hip, shoulder, extremities surgery.

It is important that the surgeon has mastered the surgical technique and that they have fully read and understood these instructions and the operating technique manual specific to the selected joint implant. The operating technique manual is provided together with the material supplied to the establishment. Alternatively, it can be consulted online or downloaded from www.evolutis-group.com upon the surgeon's request to Evolutis to obtain a personal access code for the download website. The instruments are not for single use. There is no known adverse effect when the instruments are used according to these instructions and the operating technique.

Only use cleaning agents and/or disinfectants specially formulated for cleaning or disinfecting reusable instruments recommended by Evolutis, i.e. Class I cleaning and disinfecting agents (definition of medical device according to Article 2 of Regulation 2017/745) intended for the prevention of potential infections.

The instructions provided by the detergent manufacturer regarding concentrations and temperatures must be followed.

2. LIMITATIONS OF TREATMENT

These instruments have been designed for implantation or explantation of an implant manufactured by EVOLUTIS, see the indications and contraindications in the surgical technique of the implant.

3. INITIAL TREATMENT AT THE POINT OF USE

Some instruments must be used with a motor, in this case the connection system must be a chuck or an AO tip or a Harris tip as appropriate. It is the responsibility of the operating theatre personnel to check, before the operation, the proper coupling of the instruments with the operating room equipment available. Some instruments are used to take measurements (in millimetres or degrees); the operator must ensure that the measurement is legible. Any damage to the equipment may result in malfunction. The measuring accuracy for the lengths is ± 1 mm, and of $\pm 2^\circ$ for the angles unless otherwise stated on the instruments.

The instruments should only be used for the function assigned to them in the operating techniques.

In the event of a change in the performance of the instrument, it is imperative to inform the manufacturer and quarantine the instrument.

At the end of the operation, the surgeon should pay special attention to ensure that an instrument has not been left in situ.

Instruments must not be modified.

4. MANUAL PRECLEANING

Immerse instruments completely in a Neodisher® Mediclean solution (class I medical device) of between 0.5% and 2% at 40°C maximum, and leave to soak for at least 15 minutes.

Articulated devices must be in open position, instruments with removable parts must be dismantled.

Use a soft nylon brush to gently brush the device until all visible dirt is removed for at least 3 minutes. For holes and cavities, use a swab. Activate mobile parts 10 times.

Rinse under tap water for at least 3 minutes.

5. AUTOMATIC CLEANING (1)

Place instruments in the basket of the washer-disinfector avoiding any contact between them, and if possible, orient the holes downwards.

Step	Water quality Detergent concentration range	Temperature range	Duration range
Prewash	Tap water	Cold	5 minutes
Wash	Neodisher® Mediclean type between 0.2% and 1% (class I medical device)	Between 40°C and 60°C	10 minutes
Rinse	Demineralized water	Ambiant temperature	5 minutes
Thermal disinfection	Purified water	93°C	10 minutes

(1) Treatment in a washer-disinfector according to its instructions and in accordance with ISO 15883

6. VISUAL INSPECTION

All instruments should be visually inspected to ensure that there is no contamination or damage to the instrument.

Any instrument showing signs of wear, including but not limited to: impact marks, lack of material, corrosion, any damaged surfaces, splinters of material likely to compromise its operation must be replaced.

The markings on the instruments used to measure anatomical dimensions must be legible. Inform your EVOLUTIS representative if the scales or other markings are not legible (e.g. batch number).

Drills, pins, rasps and other cutting devices should be examined carefully after treatment with alkaline detergents to ensure that their cutting edges are functional. Cutting edges should be free of notches and have a continuous edge.

For articulated instruments: check that joint movements are regular, without excessive « play ». Locking mechanisms should be checked to ensure that they work properly.

7. MAINTENANCE, INSPECTIONS AND TESTS

Instruments are subject to normal wear and tear during operations. When the operating theater staff considers that an instrument is no longer usable, they must decontaminate it, clean it and send it back to EVOLUTIS with any proof of decontamination.

8. STERILISATION

The instruments provided by Evolutis are non-sterile.

Before use, theater staff must ensure that the instruments have been sterilized. Instruments must be sterilized by trained and competent personnel.

Containers do not provide a sterile barrier: sterilization packaging must be used to maintain sterility. EVOLUTIS instruments have been validated for the following steam autoclave sterilization cycle:

- 1- Vacuuming of the sterilizer chamber 2 successive times separated by a steam injection.
- 2 - Carrying out of 5 successive pulsations.
- 3 - Steam injection.
- 4 - Sterilization at 134°C for 18 minutes, the beginning of the sterilization time countdown is indicated by a temperature probe located at the coldest zone of the autoclave.
- 5 - Steam draining.
- 6 - Vacuum drying for 30 minutes

It is essential that the sterilization method is validated for the particular type of sterilization equipment and for each load/product configuration according to ISO 17665-1.

9. PACKAGING

Clean instruments must be packaged appropriately. Special attention must be paid to cutting instruments, which must be protected.

10. STORAGE

Instruments must be stored in a clean, dry room at a temperature of 5 to 40°C (41 to 104°F).

11. BENEFITS AND PERFORMANCES

Indirect clinical benefit:

The surgical instruments allow the surgeon to properly perform the orthopaedic replacement procedure with compatible implants.

Performances claims:

- Durability of the instruments under normal conditions of use (wear resistance, resistance to mechanical stress during handling by the user, low breakage occurrence);
- Geometric and dimensional adequacy of the instruments to the implant ranges: the instruments are designed to fit the devices of which they allow implantation.

Safety claims:

Very low incidence of serious adverse events related to the use of orthopaedic surgical instruments.

12. ADDITIONAL INFORMATION

In the event of a special or unusual risk, the instruments must be disposed of according to the legislation in force.

The various materials used for the instruments comply with the standards in force.

Metallic materials are radiopaque, the plastic materials are radiotransparent.

In the event of a severe incident, the user must notify the manufacturer qualite@evolutis-group.com and the competent authority of the Member State in which the user is established.

13. MANUFACTURER CONTACT

For further information please contact EVOLUTIS

Tel 33 (0) 4 77 60 79 99 - mail qualite@evolutis-group.com or its representative.

14. COMPOSITION

- tainless steel compliant with NF EN 10088-3, NFS 94-090, NF EN 10088-2 or NF EN 10028-7 or ASTM F899, NF EN 10216-5 or NF EN 10151
- Stainless steel compliant with 5832-1 or 5832-9
- Cobalt chrome compliant with 5832-4
- Titanium compliant with 5832-3
- Polymers

The instructions provided above have been validated by the manufacturer of the medical device for preparing a medical device for reuse. It remains the responsibility of the operator to ensure that the treatment, as it is actually applied by implementing the equipment, materials and personnel of the treatment facility, achieves the desired result.

This requires verification and/or validation and routine inspection of the process.

These instructions have been drawn up in accordance with ISO 17664-1 and ISO 17664-2.

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