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GB - IMPLANTS FOR KNEE JOINT
INSTRUCTIONS FOR USE
VERSION 2025-04

For the attention of the surgeon and operation theatre staff

1. GENERAL

It is important that the surgeon becomes familiar with the surgical technique and has read these instructions and the specific surgical technique for the selected joint implant. The surgical technique is available when the equipment is delivered, or can be consulted online or downloaded from www.evolutis-group.com. As with implant placement and positioning, patient selection is an important factor. Overweight or inappropriate functional requirements can create exceptional stresses and reduce implant life. Warnings and instructions must be strictly observed.

It is specified that all implants marketed by EVOLUTIS are supplied sterile (sterilization by gamma radiation).

2. INDICATIONS - CONTRA-INDICATIONS – COMPLICATIONS

Total knee arthroplasties are indicated for the treatment of symptomatic pain and/or functional limitations of the knee in patients whose skeleton is mature and only when pain killer medication and correctly followed conservative treatment have failed. For the patient, his anatomy and the structure of his articulation will need to be adapted to receive the selected implant(s).

Rolflex TONIC is intended to restore function to the knee joint as part of a total knee arthroplasty.

2.1 Indications

The indications for ROLFLEX TONIC are:

- Degenerative non inflammatory knee disease (osteoarthritis, osteonecrosis, post-traumatic arthritis).
- Inflammatory knee disease (rheumatoid arthritis).
- Revision after failure of a conservative surgery, a ligamentoplasty, a menisectomy, an osteosynthesis, or a partial or total arthroplasty.

2.2 Limits of indication

In primary surgery of the knee joint, and even more in revision or in tumoral surgery, an important deviation of the anatomic axis of the lower limb (genu-varum or genu-valgum superior to 15° of deviation), or the quality of the bone stock and the bone defects due to the removal of any previously implanted material, can compromise the primary fixation of the implantable device and thus limit its clinical results. Depending on the degree of axial deviation, or on the location and the extension of the bone defect, a knee implant with complementary fixation means or bone loss substituting elements such as intramedullary keels or wedges and augments, and/or with a more constraint design can be considered.

2.3 Contra-indications

The arthroplasty of the knee can be contra-indicated in cases of: local or systemic infection, mental deficiency, neuromuscular afflictions, vascular or neurological afflictions, patients subject to alcohol or psychotropic problems, excessive medication, excessive functional situations (such as carrying out a sport with prevalent risk of fall or a situation in which the forces on the implant would exceed the mechanical resistance of the material), overweight, insufficient bone stock or weak demineralized bone impeding a good fixation, or severe extra articular deformation.

Implantation of prosthesis is not recommended in pregnant women.

2.4 Complications

Implanting a total knee prosthesis can lead to undesirable effects such as: infection, aseptic loosening, lysis/resorption, instability, pain, peri-prosthetic fracture, patella problems (excluding fracture and loosening), dislocation/subluxation, wear, progression of disease, implant breakage/dissociation, malalignment/malposition/incorrect sizing of the implant, arthrofibrosis, stiffness, CRPS (Complex Regional Pain Syndrome), metal related pathology including allergic reaction to prosthetic materials, femoral notching, soft-tissue problems, migration of the implant, heterotopic ossification, wound problems, limb length discrepancy and medical complications such as Deep Vein Thrombosis (DVT) or Pulmonary Embolism (PE).

2.5 Benefits and performances

Clinical benefits to the patients are pain relief, mobility restoration and function restoration.

Clinical performance and safety of the device are long term survival of the implant in situ with preservation of the function over time, and low frequency of complications.

3. PRECAUTIONS

3.1 Preoperative step

Before surgery, the surgeon must interview the patient to assess his or her physical and mental state with regard to the type of treatment planned. The surgeon must inform the patient of the implications of knee arthroplasty in everyday life with a knee prosthesis, of the pre-, intra- and post-operative steps and instructions, and of possible general and specific complications. In particular, patients should be informed of factors that could limit the performance and alter the stability or resistance of the implant(s) (e.g. excess weight, excessive functional demands), in order to delay the onset of complications that could be linked to these factors, and that they should immediately inform their surgeon of any change in the device's performance. The need to follow the surgeon's prescriptions must be clearly understood by the patient. The surgeon must ensure that the patient is able to understand and follow these instructions.

At the time of surgery, an appropriate set of implants must be available, together with the dedicated instrumentation for placing the implant(s).

3.2 Manipulation

For implants sterilized by gamma radiation, the integrity of the packaging, which guarantees the sterility of the contents, must be checked.

The expiration date of each individual package should also be checked before opening. If the packaging is damaged, do not use or re-sterilize, and return the implant in its defective packaging to EVOLUTIS or its representative. Unused components must not be re-sterilized.

Implants must be handled with extreme care to avoid scratching or damage. They should be stored in their original sealed packaging, between 5 and 40°C. If protective elements are supplied, these should be kept in place until the implants are used.

Implants are designed for single use only and must never be re-implanted.

3.3 Surgical technique

The surgeon must be fully conversant with such surgical techniques.

Careful pre-operative planning on the basis of radiographic findings should be carried out routinely.

Careful radiographic templates for each implant in a variety of magnification ratio are available from EVOLUTIS.

3.4 Fixation of implants

Coated implants must not be implanted with orthopaedic bone cement. When cement is used, it is essential that the cement manufacturer's "Instructions for Use" be followed carefully. Handling of the bone cement can influence the efficacy of the implant fixation. Care must be taken that the correct and appropriate size implant is used in conjunction with the correct instrumentation and the trial components where these are provided.

Definitive implants and trial components manufactured by Evolutis must not be used in conjunction with those of another manufacturer as components may not be compatible. The compatibility rule of the implants is described in the surgical technique of the implants through indications of sizes and match of the articular surfaces according to the type of stabilization and to the type of tibial baseplate used (fixed or mobile).

It is important to identify clearly the implant before selection and opening by a complete reading of its labeling. The labelling includes details of: stabilisation, method of fixation, side of implant, type of tibial tray, size and height.

3.5 Post operative follow up

At the end of hospitalization, the patient should be given warnings and prescriptions, preferably in writing, for complementary treatments, rehabilitation exercises and activity restrictions if necessary.

Regular follow-up with periodic check-ups is essential, and must be rigorously enforced. In all cases, tried-and-tested post-operative follow-up techniques should be followed.

In the case of MRI or CT scans, ask patients to systematically inform us that they have undergone an operation.

In the event of a serious incident, the user must notify the manufacturer and the competent authority.

3.6 Removal of implant(s)

In case of removal of the implant, the surgeon must use accepted practices with the right instrumentation.

4. WARNINGS

4.1 Re-use of an implant

The implants are intended for single use only; re-use of any implant is strictly prohibited. Should the implant be re-used its fixation and wear characteristics will be altered and the patient security put at risk.

In case of removal of an implant, this implant should be eliminated in compliance to legislation in force. In case of complete unpacking of an implant by mistake, this implant should be eliminated in compliance to legislation in force.

4.2 Imaging exam of a patient with an Evolutis prosthesis

Non-clinical tests have shown that this prosthesis is compatible with MRI under certain conditions. A patient fitted with this implant can be examined safely if the MRI examination complies with the following conditions:

- Static magnetic field equal to or less than 3 T,
- Maximum spatial magnetic field gradient of 19 T/m,
- Maximum specific absorption rate (SAR), averaged over the whole body, reported by the magnetic resonance (MR) system of 2 W/kg (normal operating mode).

Under the examination conditions defined above, the prosthesis is not expected to produce a clinically significant increase in temperature after 15 minutes of continuous scanning.

The appearance of an artifact caused by this device remains possible.

For this reason, patients wearing a prosthesis should inform the nursing staff of the presence and location of their prosthesis prior to any MRI examination.

4.3 Compatibility of implants and restrictions for use

EVOLUTIS trial implants and definitive implants must not be used with components from another manufacturer, as the components may not be compatible. The rule of compatibility between prosthetic components is mentioned in the implant operating technique, with indications of the size and compatibility of articular surfaces depending on the type of stabilisation (postero-stabilised (PS) or ultra-congruent (UC)) and the type of tibial plateau used (fixed or rotational). It is important to identify the implant selected by reading its labelling before opening and serving it. The labelling includes details of: stabilisation, method of fixation, side of implant, type of tibial tray, size and height.

4.4 Restrictions and special recommendations for Cobalt-Chromium (CoCr) components

The use of CoCr components should be avoided in:

- Patients with renal insufficiency
- Patients with cardiac insufficiency
- Patients of childbearing age.

For the other patients who will be implanted with a CoCr component, in case of suspected or proven failure of their prosthesis and whatever the component of the prosthesis involved in this failure, a medical follow-up including a cobalt blood test is recommended, knowing that cobalt is considered a CMR substance.

5. MATERIALS

The different materials used in knee implants are:

- Cobalt-Chromium: ISO standard 5832-4 or 5832-12

ISO	Cr	Mo	Ni	Fe	C	Mn	Si	W
	27 to 30	5 to 7	≤ 0.5	≤ 0.75	≤ 0.35	≤ 1.0	≤ 1.0	≤ 0.2
5832-4	P	S	N	Al	Ti	B	Co	
	≤ 0.02	≤ 0.01	≤ 0.25	≤ 0.1	≤ 0.1	≤ 0.01	Balance	

ISO	C	Cr	Ni	Mo	Mn	N	Si	Fe	Balance
5832-12	0.15 to 0.35	26 to 30	≤ 1.0	5 to 7	≤ 1.0	≤ 0.25	≤ 1.0	≤ 0.75	Co

- Polyethylene PE-UHMW (C₂H₄)_n: ISO standard 5834-1 and 2

ISO	ash	titanium	calcium	chlorine	aluminium
5834-1	≤ 125	≤ 40	≤ 5	≤ 30	≤ 20
5834-2					

- Stainless steel: ISO standard 5832-1

ISO	C	Cr	Ni	Mo	Mn	N	Si	P	S	Cu	Co	Balance
5832-1	≤ 0,03	17 to 19	13 to 15	2,25 to 3,0	≤ 2	≤ 0,1	≤ 0,75	≤ 0,025	≤ 0,01	≤ 0,50	≤ 0,10	Fe

- Pure titanium (coating): ISO standard 13179-1 and ASTM F1580

ISO	O	Fe	C	H	N	Balance
13179-1	≤ 1%	≤ 0,6%	≤ 0,1%	≤ 0,3%	≤ 0,50%	Ti

ASTM	O	Fe	C	H	N	Balance
F1580	≤ 0,40%	≤ 0,5%	≤ 0,08%	≤ 0,05%	≤ 0,05%	Ti

- T40:
- Condyles: Thickness 750µm +/-200µm in accordance with ISO 13179-1 and ASTM F1580.
 - Baseplates: Thickness from 320µm ±100µm to 140µm ±80µm in accordance with ISO 13179-1: 2021 and ASTM F1580-18.

The compositions are given in mass fraction.

The material used for each individual device is specified on the packaging of each device.

The device:

- does not incorporate a substance which, if used separately, may be considered to be a medicinal product, including a medicinal product derived from human blood or human plasma.
- is not manufactured utilizing tissues or cells of human or animal origin, or their derivatives.
- is not composed of substances or combinations of substances that are intended to be introduced into the human body and that are absorbed by or locally dispersed in the human body.
- does not contain nanomaterials.

6. IMPLANT CARD

6.1 To the attention of medical staff

An implant card is provided in the prosthesis package. You should write the patient's name and the date of the operation and attach the patient label to the space provided. This card must be given to the patient along with the information specified below.

6.2 To the attention of the patient

An implant card has been provided to you by the nursing staff following the implant of an EVOLUTIS prosthesis. You should keep this card in a safe place. It is recommended to present it during any follow-up medical visits or during an examination such as MRI, CT scan or ultrasound.

Here you will find information on the name of the implanted product, the date of operation, the manufacturer's contact details and a web link to further information on the product.

7. FURTHER INFORMATION

The summary of safety and clinical performance of device shall be made available to the public via Eudamed (<https://ec.europa.eu/tools/eudamed>). Pending the availability of the module on Eudamed, it is available on the website www.evolutis-group.com.

For any further information please contact EVOLUTIS or its local representative.

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