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GB - ORTHOPAEDIC IMPLANTS FOR SHOULDER ARTICULATION
INSTRUCTIONS FOR USE
VERSION 2024-10

For the attention of the surgeon and operation theatre staff

1. GENERAL

The operating surgeon needs to be trained and familiar with the surgical technique, and must have read and understood the IFU (Instructions For Use) leaflet and the surgical technique manual provided by EVOLUTIS with the shipment of the material to the hospital, or consulted on-line or downloaded from the www.evolutis-group.com site once the surgeon has requested and received a personal login code.

Equally to the surgical procedure and to the positioning of the implants, the selection of the patient is important. An overweight or unsuitable functional requirement can create exceptional constraints and limit the lifespan of the arthroplasty. The warnings and the instructions must be strictly observed. Precautions and instructions should be rigorously followed. It is specified that all implants marketed by EVOLUTIS are supplied sterile (sterilization by gamma radiation).

2. INDICATIONS - CONTRA-INDICATIONS – COMPLICATIONS

Hemi and total shoulder arthroplasty are indicated for the treatment of symptomatic pain and/or functional problems of the shoulder in patients whose skeleton is mature and only when pain killers and conservative treatment used correctly have failed. The anatomy and structure of the patient's articulation will need to be adapted to receive the selected implant(s).

UNIC range is intended to restore function to the shoulder joint as part of a shoulder arthroplasty (hemi or total).

2.1.Indications

The UNIC Reverse range is only indicated for total shoulder arthroplasties.

The UNIC Anatomic range is indicated for hemi or total shoulder arthroplasties.

The indications for hemi or total shoulder arthroplasty are:

- Degenerative non-inflammatory arthropathy associated or not with a functional deficiency of the rotator cuff: centred or off centred shoulder arthritis, avascular necrosis of the humeral head.
- Inflammatory arthropathy such as rheumatoid arthritis, or post traumatic arthritis.
- Post Traumatic degenerative arthritis.
- Functional repair of recent complex proximal humeral fractures.
- Revision of a previously failed hemi or total arthroplasty.

Regardless of the pathology, the presence of a functional or repairable rotator cuff is absolutely necessary for the indication of a hemi or total anatomic shoulder arthroplasty.

2.2.Contra-indications

Absolute contraindications to anatomic or reverse hemi or total shoulder arthroplasty are: local or systemic infection, mental deficiency, neuromuscular disease, neurological or vascular disease, patients addicted to alcohol or psychotropic drugs and patients with allergies to prosthetic materials. The presence of a deficient or irreparable rotator cuff is an absolute contraindication to hemi or total anatomic shoulder arthroplasty.

Relative contra indications of hemi or total shoulder arthroplasty are: excessive functional demands (sport with risk of falls, or with excessive functional demands beyond the limits of the prosthesis' resistance) overweight, insufficient bone stock or demineralisation of the bone which would compromise prosthetic fixation, severe humerus deformity, drug abuse, pre-existing peri articular oncological pathology.

Implantation of prosthesis is not recommended in pregnant women.

The UNIC Reverse range is only indicated for total shoulder arthroplasties.

The UNIC Anatomic range is indicated for hemi or total shoulder arthroplasties.

2.3.Complications

Implanting a prosthesis can lead to undesirable effects such as: infection, aseptic loosening, instability/dislocation, lysis/resorption, secondary deficiency of the rotator cuff, peri-prosthetic fracture, implant breakage/dissociation, glenoid erosion, pain, polyethylene wear, progression of disease, arthrofibrosis, metal related pathology, stiffness, soft-tissue problems, osseous tissue problems, migration of the

implant, malalignment/malposition/incorrect sizing of the implant, glenoid perforation, ossification, hematoma, wound problems and limb length discrepancy.

2.4.Benefits and performances

Clinical benefits to the patients are pain relief, mobility restoration and function restoration.

Clinical performance and safety of the device are long term survival of the implant in situ with preservation of the function over time and low frequency of complications.

3. PRECAUTIONS

3.1.Pre operative phase

Prior to the operation, the surgeon must interview the patient about his physical and mental status and difficulties and inform him about the product and the different aspects of the operation. During the consultation the surgeon should explain the limits of the operation as well as the mechanical or the material limitations of the chosen implant(s). Factors which could limit the performance of the implant and alter its stability (eg. Overweight and functional demands) should be explained in order to reduce and delay the onset of complications linked to these factors. The necessity of following the surgeons' instructions and prescription must be explained to the patient. At the time of operation the appropriate implants must be available.

3.2.Manipulation

For sterile implants the integrity of the packaging must be verified as it is a guarantee of sterility, as well as end of shelf life date. If the packaging is damaged do not use, do not re-sterilize. Implants which have been unpacked and not implanted should not be re-sterilized. Implants must be manipulated with caution in order to avoid damage or scratching. They should be conserved in their original sealed packaging between 5 and 40°C. If protective components are delivered with the implant, they must not be removed until implantation of the implant. Implants are designed for single use and must never be re-implanted.

3.3.Operative (Surgical) technique

The surgeon must use known and recognized surgical techniques. Careful pre-operative planning on x rays should be done systematically. Templates are available for the majority of our implants.

3.4.Fixation of implants

HA coated implants must not be implanted with cement. When cement is used, it is essential to follow the instructions and advice of the manufacturer. The cement mixing and introduction technique can influence the fixation of the implant. The appropriate implant of appropriate size should be used with the appropriate instrumentation and trials (if the latter are supplied).

3.5.Post operative follow up

Towards the end of the hospital stay patients should receive written advice and warnings about complimentary treatment, physio-therapy exercises and activity restrictions if necessary. In all cases validated post-operative follow up techniques must be used. The patient must be instructed to systematically inform that he or she has had an operation in case of subsequent MRI or scanner examination.

In the event of a serious incident, the user must notify the manufacturer and the competent authority.

3.6.Removal of implant(s)

In case of implant removal the surgeon must use proven surgical techniques and appropriate instrumentation.

4. WARNINGS

4.1.Re-use of an implant

Implants are single use and their re-use is forbidden. In case of re-use, the implant stability and its wear can be altered and the patient security compromised. If an implant is unwrapped by mistake or removed, they must be disposed of according to the legislation.

4.2.Imaging exam of a patient with an Evolutis prosthesis

Non-clinical testing has demonstrated that this prosthesis is MR Conditional. A patient with this implant can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 3T or less
- Maximum spatial gradient magnetic field of 19 T/m
- Maximum MR System reported, whole body averaged specific absorption rate (SAR) of 2 W/kg (normal operating mode)

Under the scan conditions defined above, the product is not expected to produce a clinically significant temperature rise after 15 minutes of continuous scanning.

The appearance of an artifact caused by this device remains possible.

That's why the patient with a prosthesis must inform the healthcare staff before any MRI examination of the presence and site of implantation of his prosthesis.

4.3. Compatibility of implants and restrictions for use

Implants and trials provided by Evolutis should not be used with components supplied by another manufacturer since they might not be compatible.

The compatibility of the implants is mentioned in the surgical technique by indications of diameters (head / insert and insert / cup association).

These characteristics are indicated on the labeling of the implants.

5. MATERIALS

The different materials used in shoulder implants are:

Stainless steel: ISO standard 5832-9 and ISO 5832-1

ISO	C	Cr	Ni	Mo	Mn	N	Si	P	S	Cu	Balance
5832-1	≤ 0,03	17 to 19	13 to 15	2,25 to 3,0	≤ 2	≤ 0,1	≤ 1	≤ 0,025	≤ 0,01	≤ 0,50	Fe

ISO	C	Cr	Ni	Mo	Mn	N	Si	P	S	Cu	Nb	Balance
5832-9	≤ 0,08	19,5 to 22,9	9 to 11	2 to 3	2 to 4,25	0,25 to 0,50	≤ 0,75	≤ 0,025	≤ 0,01	≤ 0,25	0,25 to 0,80	Fe

Titanium alloy Ti6Al4V: ISO standard 5832-3

ISO	C	N	Al	O	V	Fe	H	Balance
5832-3	≤ 0,08	≤ 0,05	5,5 to 6,75	≤ 0,2	3,5 to 4,5	≤ 0,30	≤ 0,015	Ti

Polyethylene PE-UHMW: ISO standard 5834-1 and 2

ISO	ash	titanium	calcium	chlorine	aluminium
5834-1	≤ 125	≤ 40	≤ 5	≤ 30	≤ 20

Pure titanium (coating): ISO standard 13179-1

ISO	C	N	O	Fe	H	Balance
13179-1	≤ 0,10	≤ 5	≤ 10	≤ 0,60	≤ 0,30	Ti

Calcium Hydroxyapatite (coating): ISO standard 13779-2, 3 and 4

Crystallinity: 45%

Ca/P ratio: between 1.63 and 1.824

CaO: less than 5%

The material used for each individual device is specified on the packaging of each device.

The device:

- does not incorporate a substance which, if used separately, may be considered to be a medicinal product, including a medicinal product derived from human blood or human plasma.
- is not manufactured utilising tissues or cells of human or animal origin, or their derivatives.
- is not composed of substances or combinations of substances that are intended to be introduced into the human body and that are absorbed by or locally dispersed in the human body.
- does not contain nanomaterials.

6. IMPLANT CARD

6.1. For the attention of the medical personnel

An implant card is provided in the prosthesis package. You should write the patient's name and the date of the operation and attach the patient label to the space provided. This card must be given to the patient along with the information specified below.

6.2. For the attention of the patient

An implant card has been provided to you by the nursing staff following the implant of an EVOLUTIS prosthesis. You should keep this card in a safe place. It is recommended to present it during any follow-up medical visits or during an examination such as MRI, CT scan or ultrasound.

Here you will find information on the name of the implanted product, the date of operation, the manufacturer's contact details and a web link to further information on the product.

7. FURTHER INFORMATION

The summary of safety and clinical performance of device shall be made available to the public via Eudamed (<https://ec.europa.eu/tools/eudamed>). Pending the availability of the module on Eudamed, it is available on the website www.evolutis-group.com.

For any further information please contact Evolutis or its local representative.

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