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INSTRUCTIONS FOR USE

VERSION 2024-05

For the attention of the surgeon and operation theatre staff

1. GENERAL

It is essential that the surgeon be fully conversant with the prescribed surgical technique and that due consideration be given to the information provided herewith. The specific surgical technique is delivered together with the implants and instruments to the care center before the surgery, or may be consulted online or downloaded from the www.evolutis-group.com website. The surgeon will receive a personal login to access the website once he has registered on the website. Improper patient selection, with particular regard to patient weight and functional demands, as well as placement and positioning of the implant may result in unusual stress conditions and subsequent reduction in their functional life. Warnings and instructions should be strictly observed. It is specified that all implants marketed by EVOLUTIS are supplied sterile (sterilisation by gamma radiation).

2. INDICATIONS - CONTRA-INDICATIONS - COMPLICATIONS

Hemi and total hip arthroplasties are indicated for the treatment of symptomatic pain and/or functional problems of the hip in patients whose skeleton is mature and only when pain killer medication and correctly followed conservative treatment have failed. For the patient, his anatomy and the structure of his articulation will need to be adapted to receive the selected implant(s).

2.1. Indications

Indications for hip arthroplasty are:

- Degenerative non inflammatory hip disease (coxarthrosis).
- Inflammatory hip disease (rheumatoid arthritis, spondylitis).
- Metabolic hip disease (chondrocalcinosis).
- Post Traumatic degenerative arthritis.
- Avascular necrosis.
- Congenital dysplasia of the hip.
- Functional repair of a recent trauma (fracture, dislocation)
- Revision of a failed conservative surgery, osteosynthesis, partial or total arthroplasty.
- Tumoral surgery when affecting the hip joint.

Depending on the bone condition, the indication and the approach, the use of a short, standard or long stem should be considered.

The shortened stems (STEMSYS MI) are specially designed for anterior approaches. They are not recommended for hemi-prosthesis or revision indications.

Long stems are especially dedicated to revision surgery. Stems with distal locking are recommended for advanced revision surgery or for approaches including an extended trochanterotomy or femorotomy.

The trochanteric grip plate is used in combination with the PRIUS stem and cortical screws Ø 4.5mm (see surgical technique) in cases of:

- Fracture of the trochanteric region during the placement of the PRIUS stem.
- Trochanterotomy or femorotomy necessary for the placement of the PRIUS stem or the removal of the previous femoral stem (fixation of the trochanter or the femoral flap).
- Treatment of a greater trochanter non-union secondary to the placement of a PRIUS stem.

2.2. Contra-indications

The arthroplasty of the hip can be contra-indicated in cases of: local or systemic infection, mental deficiency, neuromuscular disease, neurologic or vascular problems, patients addicted to alcohol or psychotropic drugs, drug abuse, excessive functional expectations (sport with risk of falls or with excessive functional expectations beyond the limits of the mechanical resistance of the prosthesis), overweight, insufficient bone stock or demineralisation of the bone which would compromise the prosthetic fixation, or severe extra articular deformity. Implantation of prosthesis is not recommended in pregnant women.

2.3. Complications

The implantation of a hip prosthesis may entail the following complications: Haematoma, thrombosis, pulmonary embolism, cardiovascular disorders, lack of bone union, muscle, nerve, tendon or venous damage, peri-prosthetic ossifications, allergic reactions to prosthetic materials, tissue reaction to metal or plastic

particles (metallosis-osteolysis), systemic reaction to the release of metal ions (thyroid, neurological, cardiac problems...), migration of the implant, disassembly of prosthetic modules, prosthesis-related pain, bone fracture, fracture of a prosthetic element, bone wear, prosthetic wear, prosthetic noise (snapping, squeaking), uneven limb length, dislocation, loosening of the prosthesis, infection, acute renal failure, stroke.

2.4. Benefits and performance

The intended clinical benefits and performances of these hip joint implants are:

- Pain relief,
- Improved joint mobility and functionality, contributing to a better quality of life for the patient.
- Minimum acceptable prosthesis lifetime compared to other validated prosthetic solutions (expressed as a survival rate, see acceptability threshold in the evaluation report of each prosthesis).

3. PRECAUTIONS

3.1. Preoperative step

The surgeon should discuss all physical and mental limitations particular to the patient and products for all aspects of the surgery at the pre, intra, and post-operative steps, including all the possible complications, with the patient before surgery. This discussion should include limitations of the procedure and the mechanical limitations of the material(s) selected for implantation. Factors that could impose limits on the performance and stability of the implant(s) (ex: activity level and patient weight) should be described to the patient so as to maximise the potential for long-term freedom from complications, and the patient must be instructed that any change in the performance of the implant(s) needs to be immediately reported to the surgeon for proper action. The necessity to follow the surgeon's instructions regarding these factors following the operation must be clearly understood by the patient. The surgeon must make sure that the instructions have been clearly understood by the patient. An adequate inventory of sterile implants sizes and the dedicated instrumentation set should be on hand at the time of surgery.

3.2. Handling of the implantable device(s)

For implants that are provided sterile through Gamma irradiation, the integrity of the packaging should be checked to ensure that the sterility of the contents is not compromised. The expiry date of each individual packaging should also be controlled before opening. If the packaging is damaged do not use or re-sterilise, and return the implant together with its defective packaging to Evolutis or to its local representative. Unused components should not be re-sterilized.

Extreme care must be taken while handling the implant so that it is not scratched or damaged in any way.

Implants should be stored in their original sealed packaging and in a dry and tempered room (5-40°C). Where protective covers are provided these should be left in place until immediately prior to use.

Implants are designed for single use and must never be re-implanted.

3.3. Surgical technique

The surgeon must be fully conversant with such surgical techniques.

Careful pre-operative planning on the basis of radiographic findings should be carried out routinely. Radiographic templates are available for the majority of implants.

3.4. Implant fixation

Coated implants must not be implanted with cement.

When cement is used, it is essential that the cement manufacturer's "Instructions for Use" be followed carefully. Handling of the bone cement can influence the efficacy of the implant fixation.

Care must be taken that the correct and appropriate size implant is used in conjunction with the correct instrumentation and the trial components where these are provided.

The implantation of components assembled with a Morse taper connection require extra care from the surgeon especially if one of the components is made of ceramic material. The Morse taper on the femoral stem side must be thoroughly cleaned and dried before assembly of the components. If the component is made of ceramic, turn it clock-wise before impaction. Any impaction of a ceramic component must be made through a plastic instrument. Never impact a ceramic component directly with a metal tool.

3.5. Post operative follow up

Towards the end of the hospital stay patients should receive written advice and warnings about complimentary treatment, physio-therapy exercises and activity restrictions if necessary. In all cases validated post-operative follow up techniques must be used. The patient must be instructed to systematically inform that he or she has had an operation in case of subsequent MRI or scanner examination.

In the event of a serious incident, the user must notify the manufacturer and the competent authority.

3.6. Removal of implant(s)

In case of removal of the implant, the surgeon must use accepted practices with the right instrumentation.

4. WARNINGS

4.1. Re-use of an implant

Implants are single use and their re-use is forbidden.

In case of re-use, the performance of the implant would be compromised or wear would be accelerated.

If an implant is unwrapped by mistake or removed, they must be disposed of according to the legislation.

4.2. Imaging exam of a patient with an Evolutis prosthesis

Non-clinical testing has demonstrated that this prosthesis is MR Conditional. A patient with this implant can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 3T or less
- Maximum spatial gradient magnetic field of 19 T/m
- Maximum MR System reported, whole body averaged specific absorption rate (SAR) of 2 W/kg (normal operating mode)

Under the scan conditions defined above, the product is not expected to produce a clinically significant temperature rise after 15 minutes of continuous scanning.

The appearance of an artifact caused by this device remains possible.

That's why the patient with a prosthesis must inform the healthcare staff before any MRI examination of the presence and site of implantation of his prosthesis.

4.3. Combination of implants and restrictions for use

4.3.1. Combination

By default, any femoral head (whatever the material, diameter and neck length) manufactured or distributed by EVOLUTIS is compatible with any femoral stem and any acetabular component manufactured or distributed by EVOLUTIS according to the rules of adaptation and the limits of use described in the documentation, in the IFU and on the packaging labels.

The combination of any femoral head (whatever the material, diameter and neck length) manufactured or distributed by another entity with any femoral stem manufactured or distributed by EVOLUTIS will only be acquired when EVOLUTIS has validated this combination by examining the technical file of the femoral head concerned and/or by carrying out compatibility tests.

The combination of any acetabular component manufactured or distributed by another entity with any femoral component (stem and associated head) manufactured, distributed or validated by EVOLUTIS will only be acquired when EVOLUTIS has validated this compatibility by examining the technical file of the acetabular component concerned and/or by carrying out compatibility tests.

Without validation by EVOLUTIS of the combination between the different components of a total or partial hip prosthesis including an EVOLUTIS stem or head:

- The combination of a stem manufactured or distributed by EVOLUTIS with a femoral head manufactured or distributed by another entity and
- The combination of a femoral component associating stem and head manufactured, distributed or validated by EVOLUTIS with an acetabular component manufactured or distributed by another entity are prohibited and under the sole responsibility of the prescriber and/or manufacturer or distributor concerned.

The combination of the implants is described in the surgical technique with indications for the Morse taper (combination of head/stem).

4.3.2. Restrictions

There may be restrictions on the use of certain implant sizes or models, particularly depending on their maximum mechanical strength. There are four restrictions on use:

- Weight limit to 50kg and not to be used with heads $\geq +7\text{mm}$;
- Weight limit to 70kg;
- Weight limit to 75kg;
- Weight limit to 75kg and do not use with heads $\geq +7\text{mm}$.

These restrictions shall be specified on the labels of the implants concerned.

The combination of a metal head with a conventional polyethylene insert is only recommended for patients over 70 years of age, or for patients under 70 years of age with a low life expectancy.

EVOLUTIS trial implants and definitive implants must not be used with elements from any other manufacturer, as the components may not be compatible.

These characteristics are indicated on the implant labels.

Do not combine a metal head with a ceramic insert or a ceramic head with a metal insert.

4.4. Restrictions and specific recommendations for cobalt chrome (CoCr) heads

The use of CoCr heads should be avoided in:

- Patients with renal insufficiency
- Patients with cardiac insufficiency
- Patients of childbearing age
- Patients with a known allergy to metals
- Patients undergoing revision hip surgery for which a head change is planned without a stem change, regardless of the cause of the revision

The alternatives proposed for these patients according to their level of activity and/or life expectancy are:

- If the required joint diameter is strictly less than 36mm, the use of a ceramic head or a stainless steel head if the patient has no known allergy to metals, in particular nickel,
- If the joint diameter is greater than or equal to 36mm, the use of a ceramic head.

It is recommended not to use a ceramic head in cases of hip revision if the stem taper, after careful examination, shows evidence of damage (see paragraph 4.6).

For the other patients who will be implanted with a CoCr head, in case of suspected or proven failure of their prosthesis and whatever the component of the prosthesis involved in this failure, a medical follow-up including a cobalt blood test is recommended, knowing that cobalt is considered a CMR substance.

4.5. Main risk factors for complications with metal heads

The risk of complications related to the release of metal particles from the junction between a metal head and a metal stem (combined phenomena of galvanic corrosion and mechanical fretting) is increased with the combination of the following factors

- Head and stem made of different metallic materials, in particular chrome cobalt head on titanium stem
- Head diameter strictly greater than 32mm
- Head with long or extra long neck
- Head associated with a lateralised and/or varised stem
- Patient significantly overweight (obese) and/or very active

The increase in risk depends on the number of associated factors.

Known complications of metal release can be metallosis, ALVAL (Aseptic Lymphocytic Vasculitis Associated Lesion) inflammatory reaction, premature wear of the polyethylene associated with osteolysis or prosthetic loosening.

4.6. Specific warnings in case of breakage of a ceramic component

Whenever a ceramic component breaks – except if it happens in the intraoperative step– is deemed to be destructive for the metallic taper on which it is fixed. Consequently:

- In case of breakage of the ceramic femoral head when used in a ceramic-polyethylene articulation, the general guidance is to replace the femoral stem and the broken femoral head by a new femoral stem and a new ceramic head, except if – after a thorough inspection of the connection taper of the femoral stem – the taper is not damaged; in this case, the replacement of the femoral head alone may be considered. In any case, it is recommended to change the polyethylene insert.
- In case of breakage of one of the ceramic components when used in a ceramic-ceramic articulation, the general guidance is to replace the 2 ceramic articular components, the femoral stem and the acetabular cup. An exception can be made for the articular component which is not damaged if – after removal of the ceramic component and thorough inspection of the connection taper – the taper is not damaged: the replacement of the ceramic component alone can be considered.
- In case of breakage of the two ceramic components when used in a ceramic-ceramic articulation, the general guidance is to replace the 2 ceramic articular components, the femoral stem and the acetabular cup.

In all cases of breakage of a ceramic component, it is recommended to thoroughly clean the implantation site in order to remove any remaining ceramic fragment likely to induce third body complication.

5. MATERIALS

The different materials used in hip implants are:

- Stainless steel: ISO standard 5832-1 or 5832-9

ISO	C	Cr	Ni	Mo	Mn	N	Si	P	S	Cu	Other
5832-1	≤ 0,03	17 to 19	13 to 15	2,25 to 3,0	≤ 2	≤ 0,1	≤ 1	≤ 0,025	≤ 0,01	≤ 0,50	Fe

ISO	C	Cr	Ni	Mo	Mn	N	Si	P	S	Cu	Nb	Other
5832-9	≤ 0,08	19,5 to 22	9 to 11	2 to 3	2 to 4,25	0,25 to 0,50	≤ 0,75	≤ 0,025	≤ 0,01	≤ 0,25	0,25 to 0,80	Fe

- Titanium alloy: ISO standard 5832-3

ISO	C	N	Al	O	V	Fe	H	Other
5832-3	≤ 0,08	≤ 0,05	5.5 to 6.75	≤ 0,2	3.5 to 4.5	≤ 0.30	≤ 0.015	Ti

- Cobalt chrome: ISO standard 5832-4 or 5832-12

ISO	C	Cr	Ni	Mo	Mn	N	Si	Fe	Other
5832-4	0,15 to 0,35	26 to 30	≤ 1,0	5 to 7	≤ 1,0	≤ 0,25	≤ 1,0	≤ 1,0	Co

ISO	C	Cr	Ni	Mo	Mn	N	Si	Fe	Other
5832-12	≤ 0,14	26 to 30	≤ 1,0	5 to 7	≤ 1,0	≤ 0,25	≤ 1,0	≤ 0,75	Co

- Ceramic: ISO standard 6474-2

ISO	ZrO2	HfO2	Impurity	Additive	Other
6474-2	10 to 30	≤ 5	≤ 0,2	≤ 10	Al2O3

- Pure titanium (coating): ISO standard 13179-1

ISO	C	N	O	Fe	H	Other
13179-1	≤ 0,10	≤ 5	≤ 10	≤ 0,60	≤ 0,30	Ti

- Calcium Hydroxyapatite (coating): ISO standard 13779-2, 3 and 4

Ca10 (PO4)6 (OH)2 in compliance to ISO 13779-2, -3 and -4.

Crystallinity: 45%

Ca/P ratio: between 1.63 and 1.824

CaO: less than 5%

The compositions are given in mass fraction.

The material used for each individual device is specified on the packaging of each device.

6. IMPLANT CARD

6.1. For the attention of the medical personnel

An implant card is provided in the prosthesis package. You should write the patient's name and the date of the operation and attach the patient label to the space provided. This card must be given to the patient along with the information specified below.

6.2. For the attention of the patient

An implant card has been provided to you by the nursing staff following the implant of an EVOLUTIS prosthesis. You should keep this card in a safe place. It is recommended to present it during any follow-up medical visits or during an examination such as MRI, CT scan or ultrasound.

Here you will find information on the name of the implanted product, the date of operation, the manufacturer's contact details and a web link to further information on the product.

7. FURTHER INFORMATION

The summary of safety and clinical performance of device shall be made available to the public via Eudamed (<https://ec.europa.eu/tools/eudamed>). Pending the availability of the module on Eudamed, it is available on the website www.evolutis-group.com.

For any further information please contact Evolutis or its local representative.

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